



CACNA1B Antibody, FITC conjugated

Product Code	CSB-PA004398LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00975
Immunogen	Recombinant Human Voltage-dependent N-type calcium channel subunit alpha-1B protein (2132-2339AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent N-type calcium channel subunit alpha-1B (Brain calcium channel III) (BIII) (Calcium channel, L type, alpha-1 polypeptide isoform 5) (Voltage-gated calcium channel subunit alpha Cav2.2), CACNA1B, CACH5 CACNL1A5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CACNA1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.