



CD81 Antibody, Biotin conjugated

Product Code	CSB-PA004960ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60033
Immunogen	Peptide sequence from Human CD81 antigen protein (116-134AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD81 antigen (26 kDa cell surface protein TAPA-1) (Target of the antiproliferative antibody 1) (Tetraspanin-28) (Tspan-28) (CD antigen CD81), CD81, TAPA1 TSPAN28
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD81
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.