




DVL1P1 Antibody

Product Code	CSB-PA007285OA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54792
Immunogen	Peptide sequence from Human Putative segment polarity protein dishevelled homolog DVL1P1 protein (152-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative segment polarity protein dishevelled homolog DVL1P1 (DSH homolog 1-like) (Segment polarity protein dishevelled homolog DVL-1-like) (Dishevelled-1-like), DVL1P1, DVL DVL1 DVL1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	DVL1P1

Image



Western Blot

Positive WB detected in: PC-3 whole cell lysate

All lanes: DVL1P1 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 74, 73 kDa

Observed band size: 74 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.