



DVL1P1 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA007285OC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P54792 |
| Immunogen | Peptide sequence from Human Putative segment polarity protein dishevelled homolog DVL1P1 protein (152-173AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Putative segment polarity protein dishevelled homolog DVL1P1 (DSH homolog 1-like) (Segment polarity protein dishevelled homolog DVL-1-like) (Dishevelled-1-like), DVL1P1, DVL DVL1 DVL1L1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | DVL1P1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |