



DVL1P1 Antibody, FITC conjugated

Product Code	CSB-PA007285OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54792
Immunogen	Peptide sequence from Human Putative segment polarity protein dishevelled homolog DVL1P1 protein (152-173AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative segment polarity protein dishevelled homolog DVL1P1 (DSH homolog 1-like) (Segment polarity protein dishevelled homolog DVL-1-like) (Dishevelled-1-like), DVL1P1, DVL DVL1 DVL1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	DVL1P1