



EIF4G2 Antibody, HRP conjugated

Product Code	CSB-PA007569LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78344
Immunogen	Recombinant Human Eukaryotic translation initiation factor 4 gamma 2 protein (706-907AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eukaryotic translation initiation factor 4 gamma 2 (eIF-4-gamma 2) (eIF-4G 2) (eIF4G 2) (Death-associated protein 5) (DAP-5) (p97), EIF4G2, DAP5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	EIF4G2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.