



# ERC1 Antibody, HRP conjugated

Product Code	CSB-PA007767LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUD2
Immunogen	Recombinant Human ELKS/Rab6-interacting/CAST family member 1 protein (1-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ELKS/Rab6-interacting/CAST family member 1 (ERC-1) (Rab6-interacting protein 2), ERC1, ELKS KIAA1081 RAB6IP2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ERC1