



ERC1 Antibody, HRP conjugated

Product Code	CSB-PA007767LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUD2
Immunogen	Recombinant Human ELKS/Rab6-interacting/CAST family member 1 protein (1-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ELKS/Rab6-interacting/CAST family member 1 (ERC-1) (Rab6-interacting protein 2), ERC1, ELKS KIAA1081 RAB6IP2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ERC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.