



GPR37 Antibody, Biotin conjugated

Product Code	CSB-PA009810ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15354
Immunogen	Peptide sequence from Human Prosaposin receptor GPR37 protein (114-133AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prosaposin receptor GPR37 (Endothelin B receptor-like protein 1) (ETBR-LP-1) (G-protein coupled receptor 37) (Parkin-associated endothelin receptor-like receptor) (PAELR), GPR37
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR37
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.