



GPR55 Antibody, Biotin conjugated

Product Code	CSB-PA009820OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2T6
Immunogen	Peptide sequence from Human G-protein coupled receptor 55 protein (203-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled receptor 55, GPR55
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR55
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.