




# HSPA12A Antibody

<b>Product Code</b>	CSB-PA010817LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O43301
<b>Immunogen</b>	Recombinant Human Heat shock 70 kDa protein 12A protein (2-231AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Heat shock 70 kDa protein 12A, HSPA12A, KIAA0417
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	HSPA12A

## Image



### Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate, 293 whole cell lysate, HepG2 whole cell lysate

All lanes: HSPA12A antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 75 kDa

Observed band size: 75 kDa