




HSPA12A Antibody

Product Code	CSB-PA010817LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43301
Immunogen	Recombinant Human Heat shock 70 kDa protein 12A protein (2-231AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein 12A, HSPA12A, KIAA0417
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPA12A

Image



Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate, 293 whole cell lysate, HepG2 whole cell lysate

All lanes: HSPA12A antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 75 kDa

Observed band size: 75 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.