



HSPA12A Antibody, HRP conjugated

Product Code	CSB-PA010817LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43301
Immunogen	Recombinant Human Heat shock 70 kDa protein 12A protein (2-231AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein 12A, HSPA12A, KIAA0417
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPA12A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.