



HSPA12A Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA010817LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O43301 |
| Immunogen | Recombinant Human Heat shock 70 kDa protein 12A protein (2-231AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Heat shock 70 kDa protein 12A, HSPA12A, KIAA0417 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | HSPA12A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |