



IL13RA1 Antibody, HRP conjugated

Product Code	CSB-PA011591LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78552
Immunogen	Recombinant Human Interleukin-13 receptor subunit alpha-1 protein (22-234AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-13 receptor subunit alpha-1 (IL-13 receptor subunit alpha-1) (IL-13R subunit alpha-1) (IL-13R-alpha-1) (IL-13RA1) (Cancer/testis antigen 19) (CT19) (CD antigen CD213a1), IL13RA1, IL13R IL13RA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL13RA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.