



KALRN Antibody, HRP conjugated

Product Code	CSB-PA011981LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60229
Immunogen	Recombinant Human Kalirin protein (2410-2661AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kalirin (EC 2.7.11.1) (Huntingtin-associated protein-interacting protein) (Protein Duo) (Serine/threonine-protein kinase with Dbl- and pleckstrin homology domain), KALRN, DUET DUO HAPIP TRAD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KALRN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.