



KALRN Antibody, HRP conjugated

Product Code	CSB-PA011981LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60229
Immunogen	Recombinant Human Kalirin protein (2410-2661AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kalirin (EC 2.7.11.1) (Huntingtin-associated protein-interacting protein) (Protein Duo) (Serine/threonine-protein kinase with Dbl- and pleckstrin homology domain), KALRN, DUET DUO HAPIP TRAD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KALRN