



LCHN Antibody, Biotin conjugated

Product Code	CSB-PA012231LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4D1U4
Immunogen	Recombinant Human Protein LCHN protein (2-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein LCHN, LCHN, KIAA1147
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LCHN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.