



LAMB4 Antibody, Biotin conjugated

Product Code	CSB-PA012734LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4D0S4
Immunogen	Recombinant Human Laminin subunit beta-4 protein (1513-1761AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Laminin subunit beta-4 (Laminin beta-1-related protein), LAMB4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LAMB4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.