



LPAR2 Antibody, HRP conjugated

Product Code	CSB-PA013048OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBW0
Immunogen	Peptide sequence from Human Lysophosphatidic acid receptor 2 protein (8-27AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysophosphatidic acid receptor 2 (LPA receptor 2) (LPA-2) (Lysophosphatidic acid receptor Edg-4), LPAR2, EDG4 LPA2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LPAR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.