



MCF2 Antibody, FITC conjugated

Product Code	CSB-PA013579LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10911
Immunogen	Recombinant Human Proto-oncogene DBL protein (705-925AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proto-oncogene DBL (Proto-oncogene MCF-2) [Cleaved into: MCF2-transforming protein; DBL-transforming protein], MCF2, DBL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MCF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.