



NPFFR2 Antibody, Biotin conjugated

Product Code	CSB-PA015983ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5X5
Immunogen	Peptide sequence from Human Neuropeptide FF receptor 2 protein (23-43AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropeptide FF receptor 2 (G-protein coupled receptor 74) (G-protein coupled receptor HLWAR77) (Neuropeptide G-protein coupled receptor), NPFFR2, GPR74 NPFF2 NPGPR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NPFFR2