



OPRK1 Antibody, HRP conjugated

Product Code	CSB-PA016359OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41145
Immunogen	Peptide sequence from Human Kappa-type opioid receptor protein (31-50AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kappa-type opioid receptor (K-OR-1) (KOR-1), OPRK1, OPRK
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OPRK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.