



OPTC Antibody, HRP conjugated

Product Code	CSB-PA016362LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBM4
Immunogen	Recombinant Human Opticin protein (20-332AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Opticin (Oculoglycan), OPTC, OPT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OPTC