



OPTC Antibody, HRP conjugated

Product Code	CSB-PA016362LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBM4
Immunogen	Recombinant Human Opticin protein (20-332AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Opticin (Oculoglycan), OPTC, OPT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OPTC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.