




PAX8 Antibody

Product Code	CSB-PA017494LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06710
Immunogen	Recombinant Human Paired box protein Pax-8 protein (237-450AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paired box protein Pax-8, PAX8
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PAX8

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, HeLa whole cell lysate, HepG2 whole cell lysate, PC-3 whole cell lysate

All lanes: PAX8 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 49, 42, 44, 35, 32 kDa

Observed band size: 49 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.