



PRKACA Antibody, HRP conjugated

Product Code	CSB-PA018688LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17612
Immunogen	Recombinant Human cAMP-dependent protein kinase catalytic subunit alpha protein (233-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cAMP-dependent protein kinase catalytic subunit alpha (PKA C-alpha) (EC 2.7.11.11), PRKACA, PKACA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKACA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.