



RRAS Antibody, FITC conjugated

Product Code	CSB-PA020513PC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10301
Immunogen	Peptide sequence from Human Ras-related protein R-Ras protein (7-23AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein R-Ras (p23), RRAS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RRAS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.