



TAC1 Antibody, HRP conjugated

Product Code	CSB-PA023062LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20366
Immunogen	Recombinant Human Protachykinin-1 protein (20-129AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protachykinin-1 (PPT) [Cleaved into: Substance P; Neurokinin A (NKA) (Neuromedin L) (Substance K); Neuropeptide K (NPK); Neuropeptide gamma; C-terminal-flanking peptide], TAC1, NKA NKNA TAC2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.