



ZFYVE9 Antibody, HRP conjugated

Product Code	CSB-PA026477LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95405
Immunogen	Recombinant Human Zinc finger FYVE domain-containing protein 9 protein (683-926AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger FYVE domain-containing protein 9 (Mothers against decapentaplegic homolog-interacting protein) (Madh-interacting protein) (Novel serine protease) (NSP) (Receptor activation anchor) (hSARA) (Smad anchor for receptor activation), ZFYVE9, MADHIP SARA SMADIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ZFYVE9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.