





## SOS2 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA22327C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07890
Immunogen	Recombinant Human Son of sevenless homolog 2 protein (187-404AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Son of sevenless homolog 2 (SOS-2), SOS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SOS2