



SOS2 Antibody, Biotin conjugated

Product Code	CSB-PA22327D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07890
Immunogen	Recombinant Human Son of sevenless homolog 2 protein (187-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Son of sevenless homolog 2 (SOS-2), SOS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SOS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.