



RHBDD3 Antibody, HRP conjugated

Product Code	CSB-PA27742B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3P4
Immunogen	Peptide sequence from Human Rhomboid domain-containing protein 3 protein (304-323AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rhomboid domain-containing protein 3, RHBDD3, C22orf3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RHBDD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.