



CASTOR2 Antibody, Biotin conjugated

Product Code	CSB-PA28239D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NHX0
Immunogen	Recombinant Human Cytosolic arginine sensor for mTORC1 subunit 2 protein (82-329AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic arginine sensor for mTORC1 subunit 2 (Cellular arginine sensor for mTORC1 protein 2) (GATS-like protein 2), CASTOR2, GATSL1 GATSL2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CASTOR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.