



LRRC3C Antibody, Biotin conjugated

Product Code	CSB-PA408988LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NJW4
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 3C protein (42-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing protein 3C, LRRC3C
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LRRC3C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.