



RAPGEF1 Antibody, HRP conjugated

Product Code	CSB-PA613430NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13905
Immunogen	Peptide sequence from Human Rap guanine nucleotide exchange factor 1 protein (22-40AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rap guanine nucleotide exchange factor 1 (CRK SH3-binding GNRP) (Guanine nucleotide-releasing factor 2) (Protein C3G), RAPGEF1, GRF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAPGEF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.