



PLA2G4E Antibody, HRP conjugated

Product Code	CSB-PA666060LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3MJ16
Immunogen	Recombinant Human Cytosolic phospholipase A2 epsilon protein (604-868AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic phospholipase A2 epsilon (cPLA2-epsilon) (EC 3.1.1.4) (Phospholipase A2 group IVE), PLA2G4E
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLA2G4E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.