



USP17L6P Antibody, Biotin conjugated

Product Code	CSB-PA757690LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6QN14
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 17-like protein 6 protein (136-398AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 17-like protein 6 (EC 3.4.19.12) (Deubiquitinating enzyme 17-like protein 6) (Ubiquitin thioesterase 17-like protein 6) (Ubiquitin-specific-processing protease 17-like protein 6), USP17L6P, USP17C USP17D USP17N
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	USP17L6P
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.