



THOC6 Antibody, Biotin conjugated

Product Code	CSB-PA769793LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86W42
Immunogen	Recombinant Human THO complex subunit 6 homolog protein (56-315AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	THO complex subunit 6 homolog (Functional spliceosome-associated protein 35) (fSAP35) (WD repeat-containing protein 58), THOC6, WDR58
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	THOC6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.