



VPS36 Antibody, FITC conjugated

Product Code	CSB-PA801240OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VN1
Immunogen	Peptide sequence from Human Vacuolar protein-sorting-associated protein 36 protein (87-105AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein-sorting-associated protein 36 (ELL-associated protein of 45 kDa) (ESCRT-II complex subunit VPS36), VPS36, C13orf9 EAP45
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	VPS36