



NDRG1 Antibody, FITC conjugated

Product Code	CSB-PA835678LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92597
Immunogen	Recombinant Human Protein NDRG1 protein (136-394AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein NDRG1 (Differentiation-related gene 1 protein) (DRG-1) (N-myc downstream-regulated gene 1 protein) (Nickel-specific induction protein Cap43) (Reducing agents and tunicamycin-responsive protein) (RTP) (Rit42), NDRG1, CAP43 DRG1 RTP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NDRG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.