




ADCY10 Antibody

Product Code	CSB-PA839393LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PN6
Immunogen	Recombinant Human Adenylate cyclase type 10 protein (897-1136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenylate cyclase type 10 (EC 4.6.1.1) (AH-related protein) (Adenylate cyclase homolog) (Germ cell soluble adenyl cyclase) (hsAC) (sAC) (Testicular soluble adenyl cyclase), ADCY10, SAC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADCY10

Image



IHC image of CSB-PA839393LA01HU diluted at 1:100 and staining in paraffin-embedded human testis tissue performed on a Leica BondTM system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.