



ADCY10 Antibody, FITC conjugated

Product Code	CSB-PA839393LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PN6
Immunogen	Recombinant Human Adenylate cyclase type 10 protein (897-1136AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenylate cyclase type 10 (EC 4.6.1.1) (AH-related protein) (Adenylate cyclase homolog) (Germ cell soluble adenylyl cyclase) (hsAC) (sAC) (Testicular soluble adenylyl cyclase), ADCY10, SAC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADCY10