



# ADCY10 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA839393LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96PN6
<b>Immunogen</b>	Recombinant Human Adenylate cyclase type 10 protein (897-1136AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Adenylate cyclase type 10 (EC 4.6.1.1) (AH-related protein) (Adenylate cyclase homolog) (Germ cell soluble adenylyl cyclase) (hsAC) (sAC) (Testicular soluble adenylyl cyclase), ADCY10, SAC
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	ADCY10
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.