




# CDK5RAP3 Antibody


<b>Product Code</b>	CSB-PA842691LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96JB5
<b>Immunogen</b>	Recombinant Human CDK5 regulatory subunit-associated protein 3 protein (42-292AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Rat
<b>Tested Applications</b>	ELISA, WB, IF; Recommended dilution: WB:1:500-1:2000, IF:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	CDK5 regulatory subunit-associated protein 3 (CDK5 activator-binding protein C53) (LXXLL/leucine-zipper-containing ARF-binding protein) (Protein HSF-27), CDK5RAP3, IC53 LZAP
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	CDK5RAP3

## Image



### Western Blot

Positive WB detected in: U87 whole cell lysate, Jurkat whole cell lysate, Rat liver tissue  
 All lanes: CDK5RAP3 antibody at 1:1000  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/50000 dilution  
 Predicted band size: 57, 47, 32, 60 kDa  
 Observed band size: 57 kDa



Immunofluorescence staining of HepG2 cells with CSB-PA842691LA01HU at 1:50, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.