



VSTM2L Antibody, Biotin conjugated

Product Code	CSB-PA842721LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96N03
Immunogen	Recombinant Human V-set and transmembrane domain-containing protein 2-like protein (25-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	V-set and transmembrane domain-containing protein 2-like protein, VSTM2L, C20orf102
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	VSTM2L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.