



DCHS1 Antibody, Biotin conjugated

Product Code	CSB-PA850320OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JQ0
Immunogen	Peptide sequence from Human Protocadherin-16 protein (2964-2981AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin-16 (Cadherin-19) (Cadherin-25) (Fibroblast cadherin-1) (Protein dachsous homolog 1), DCHS1, CDH19 CDH25 FIB1 KIAA1773 PCDH16
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DCHS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.