



PPP1R13L Antibody, HRP conjugated

Product Code	CSB-PA851965OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUF5
Immunogen	Peptide sequence from Human RelA-associated inhibitor protein (83-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	RelA-associated inhibitor (Inhibitor of ASPP protein) (Protein iASPP) (NFkB-interacting protein 1) (PPP1R13B-like protein), PPP1R13L, IASPP NKIP1 PPP1R13BL RAI
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PPP1R13L