



PPP1R13L Antibody, FITC conjugated

Product Code	CSB-PA851965OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUF5
Immunogen	Peptide sequence from Human RelA-associated inhibitor protein (83-102AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	RelA-associated inhibitor (Inhibitor of ASPP protein) (Protein iASPP) (NFkB-interacting protein 1) (PPP1R13B-like protein), PPP1R13L, IASPP NKIP1 PPP1R13BL RAI
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PPP1R13L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.