



CSRNP3 Antibody, HRP conjugated

Product Code	CSB-PA855508NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WYN3
Immunogen	Peptide sequence from Human Cysteine/serine-rich nuclear protein 3 protein (518-536AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cysteine/serine-rich nuclear protein 3 (CSRNP-3) (Protein FAM130A2) (TGF-beta-induced apoptosis protein 2) (TAIP-2), CSRNP3, FAM130A2 TAIP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CSRNP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.